

## STANDARD DRESSING FORCEPS

straight; serrated; narrow

20.0cm (8")

**SKU:** 020-004-20

**Categories:** [Open Surgery](#), [Forceps](#)



## ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Dressing Forceps</a> , <a href="#">Tissue Forceps</a>
<b>Aesculap</b>	BD031
<b>KLS Martin</b>	12-120-20
<b>AUTHORS</b>	<a href="#">STANDARD</a>
<b>Medicon</b>	06.01.20
<b>SPECIFICATIONS</b>	<a href="#">serrated</a> , <a href="#">straight</a>



## ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Dressing Forceps</a> , <a href="#">Tissue Forceps</a>
<b>Aesculap</b>	BD031
<b>KLS Martin</b>	12-120-20

## AUTHORS

[STANDARD](#)

## Medicon

06.01.20

## SPECIFICATIONS

[serrated](#), [straight](#)