

STANDARD DRESSING FORCEPS

straight; serrated; narrow

30.0cm (12")

SKU: 020-004-30
Categories: [Forceps](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dressing Forceps , Tissue Forceps
Aesculap	BD033
KLS Martin	12-120-30
AUTHORS	MODELL USA
Medicon	06.01.30
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dressing Forceps , Tissue Forceps
Aesculap	BD033
KLS Martin	12-120-30

AUTHORS

[MODELL USA](#)

Medicon

06.01.30

SPECIFICATIONS

[serrated](#), [straight](#)