

STANDARD DRESSING FORCEPS

straight; serrated

11.5cm (4.6"); width of jaw 2.8 +/-0.2

SKU: 020-002-11
 Categories: [Open Surgery](#), [Forceps](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dressing Forceps , Tissue Forceps
Aesculap	BD043
KLS Martin	12-100-11
AUTHORS	STANDARD
Medicon	06.00.11
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dressing Forceps , Tissue Forceps
Aesculap	BD043
KLS Martin	12-100-11

AUTHORS

[STANDARD](#)

Medicon

06.00.11

SPECIFICATIONS

[serrated](#), [straight](#)