

STANDARD DRESSING FORCEPS

straight; serrated; fine

14.5cm (5.8"); width of jaw 1.7 +/-0.2

SKU: 020-010-14
Categories: [Open Surgery](#), [Forceps](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dressing Forceps , Tissue Forceps
Aesculap	BD305
KLS Martin	12-126-14
AUTHORS	STANDARD
Medicon	06.02.14
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dressing Forceps , Tissue Forceps
Aesculap	BD305
KLS Martin	12-126-14

AUTHORS

[STANDARD](#)

Medicon

06.02.14

SPECIFICATIONS

[serrated](#), [straight](#)