

STANDARD DRESSING FORCEPS

straight; serrated; fine

20.0cm (8"); width of jaw 1.8 +/-0.2

SKU: 020-010-20
 Categories: [Open Surgery](#), [Forceps](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dressing Forceps , Tissue Forceps
Aesculap	BD031
AUTHORS	STANDARD
Medicon	06.02.20
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dressing Forceps , Tissue Forceps
Aesculap	BD031
AUTHORS	STANDARD
Medicon	06.02.20

SPECIFICATIONS

[serrated](#), [straight](#)