

STANDARD DRESSING FORCEPS

straight; serrated; fine

30.0cm (12"); width of jaw 2.5 +/-0.2

SKU: 020-010-30
Categories: [Open Surgery](#), [Forceps](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dressing Forceps , Tissue Forceps
AUTHORS	MODELL USA
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dressing Forceps , Tissue Forceps
AUTHORS	MODELL USA
SPECIFICATIONS	serrated , straight