

STANDARD DRESSING FORCEPS

straight; serrated

13.0cm (5.2"); width of jaw 2.7 +/- 0.2

SKU: 020-002-13
Categories: [Open Surgery](#), [Forceps](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dressing Forceps , Tissue Forceps
Aesculap	BD045
KLS Martin	12-100-13
AUTHORS	STANDARD
Medicon	06.00.13
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dressing Forceps , Tissue Forceps
Aesculap	BD045
KLS Martin	12-100-13

AUTHORS

[STANDARD](#)

Medicon

06.00.13

SPECIFICATIONS

[serrated](#), [straight](#)