

**STANDARD DRESSING FORCEPS**

straight; serrated

16.0cm (6.4"); width of jaw 3.7 +/-0.2

**SKU:** 020-002-16  
**Categories:** [Forceps](#), [Open Surgery](#)



**ZUSÄTZLICHE INFORMATIONEN**

|                       |                                                                   |
|-----------------------|-------------------------------------------------------------------|
| <b>UNIT</b>           | <a href="#">pcs.</a>                                              |
| <b>TYPE</b>           | <a href="#">Dressing Forceps</a> , <a href="#">Tissue Forceps</a> |
| <b>Aesculap</b>       | BD049                                                             |
| <b>KLS Martin</b>     | 12-100-16                                                         |
| <b>AUTHORS</b>        | <a href="#">STANDARD</a>                                          |
| <b>Medicon</b>        | 06.00.16                                                          |
| <b>SPECIFICATIONS</b> | <a href="#">serrated</a> , <a href="#">straight</a>               |



**ZUSÄTZLICHE INFORMATIONEN**

|                   |                                                                   |
|-------------------|-------------------------------------------------------------------|
| <b>UNIT</b>       | <a href="#">pcs.</a>                                              |
| <b>TYPE</b>       | <a href="#">Dressing Forceps</a> , <a href="#">Tissue Forceps</a> |
| <b>Aesculap</b>   | BD049                                                             |
| <b>KLS Martin</b> | 12-100-16                                                         |

## AUTHORS

[STANDARD](#)

## Medicon

06.00.16

## SPECIFICATIONS

[serrated](#), [straight](#)