

## STANDARD DRESSING FORCEPS

straight; serrated

18.0cm (7.2"); width of jaw 3.6 +/-0.2

**SKU:** 020-002-18

**Categories:** [Open Surgery](#), [Forceps](#)



## ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Dressing Forceps</a> , <a href="#">Tissue Forceps</a>
<b>Aesculap</b>	BD050
<b>KLS Martin</b>	12-100-18
<b>AUTHORS</b>	<a href="#">STANDARD</a>
<b>Medicon</b>	06.00.18
<b>SPECIFICATIONS</b>	<a href="#">serrated</a> , <a href="#">straight</a>



## ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Dressing Forceps</a> , <a href="#">Tissue Forceps</a>
<b>Aesculap</b>	BD050
<b>KLS Martin</b>	12-100-18

## AUTHORS

[STANDARD](#)

## Medicon

06.00.18

## SPECIFICATIONS

[serrated](#), [straight](#)