

STANDARD DRESSING FORCEPS

straight; serrated

20.0cm (8"); width of jaw 3.5 +/-0.2

SKU: 020-002-20

Categories: [Open Surgery](#), [Forceps](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dressing Forceps , Tissue Forceps
Aesculap	BD051
KLS Martin	12-100-20
AUTHORS	STANDARD
Medicon	06.00.20
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dressing Forceps , Tissue Forceps
Aesculap	BD051
KLS Martin	12-100-20

AUTHORS

[STANDARD](#)

Medicon

06.00.20

SPECIFICATIONS

[serrated](#), [straight](#)