

STANDARD DRESSING FORCEPS

straight; serrated

30.0cm (12"); width of jaw 5.0 +/-0.2

SKU: 020-002-30
Categories: [Forceps](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dressing Forceps , Tissue Forceps
Aesculap	BD053
KLS Martin	12-100-30
AUTHORS	STANDARD
Medicon	06.00.30
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dressing Forceps , Tissue Forceps
Aesculap	BD053
KLS Martin	12-100-30

AUTHORS

[STANDARD](#)

Medicon

06.00.30

SPECIFICATIONS

[serrated](#), [straight](#)