

STANDARD DRESSING FORCEPS

straight; serrated

10.5cm (4.2"); width of jaw 2.4 +/-0.2

SKU: 020-002-10
Categories: [Open Surgery](#), [Forceps](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dressing Forceps , Tissue Forceps
KLS Martin	12-100-10
AUTHORS	STANDARD
Medicon	06.00.10
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dressing Forceps , Tissue Forceps
KLS Martin	12-100-10
AUTHORS	STANDARD
Medicon	06.00.10

SPECIFICATIONS

[serrated](#), [straight](#)