

STANDARD FIBER COLD LIGHT GUIDE CABLE

without connectors

Ø 3,5 x 2300 mm Length

Tube colour grey; careful with Xenon light

STANDARD
COLD LIGHT LINE



SKU: 540-803-23

Categories: [Endoscopy](#), [Equipment](#), [Laparoscopy](#), [Gynecology](#), [GUBBINI/Gyn](#), [Urology](#), [Arthroscopy](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Light Guide Cables
TYPE	Standard (Cold Light)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Light Guide Cables
TYPE	Standard (Cold Light)