

STANDARD FIBER COLD LIGHT GUIDE CABLE

without connectors

Ø 4,9 x 3000 mm Length, Tube colour grey

careful with Xenon lights

STANDARD
COLD LIGHT LINE



SKU: 540-805-30

Categories: [Arthroscopy](#), [Endoscopy](#), [Equipment](#), [GUBBINI/Gyn](#), [Gynecology](#), [Laparoscopy](#), [Urology](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Light Guide Cables
TYPE	Standard (Cold Light)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Light Guide Cables
TYPE	Standard (Cold Light)