

## STANDARD FIBER COLD LIGHT GUIDE CABLE

without connectors

Ø 3,5 x 1800 mm Length

Tube colour grey; careful with Xenon light

**STANDARD**  
COLD LIGHT LINE



**SKU:** 540-803-18

**Categories:** [Arthroscopy](#), [Endoscopy](#), [Equipment](#), [GUBBINI/Gyn](#), [Gynecology](#), [Laparoscopy](#), [Urology](#)

### ZUSÄTZLICHE INFORMATIONEN

|              |                                       |
|--------------|---------------------------------------|
| <b>UNIT</b>  | <a href="#">pcs.</a>                  |
| <b>RANGE</b> | <a href="#">Light Guide Cables</a>    |
| <b>TYPE</b>  | <a href="#">Standard (Cold Light)</a> |

### ZUSÄTZLICHE INFORMATIONEN

|              |                                       |
|--------------|---------------------------------------|
| <b>UNIT</b>  | <a href="#">pcs.</a>                  |
| <b>RANGE</b> | <a href="#">Light Guide Cables</a>    |
| <b>TYPE</b>  | <a href="#">Standard (Cold Light)</a> |