

STANDARD SHEATH;17FR.; ROTABLE; 1 STOPCOC

for bipolar working element 582-180-00

QL connection, colour code green



SKU: 582-180-15

Categories: [Endoscopy](#), [Urology](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Mini Resectoscopy
TYPE	Sheaths

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Mini Resectoscopy
TYPE	Sheaths