

STEINMANN EXTENSION PIN

with trocar point; triangular shaft

Ø 4.0mm; 20.0cm (8")



SKU: 212-348-20

Categories: [Open Surgery](#), [Orthopaedics](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Steinmann pins
Aesculap	KX120
KLS Martin	22-562-20
AUTHORS	STEINMANN

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Steinmann pins
Aesculap	KX120
KLS Martin	22-562-20
AUTHORS	STEINMANN