

STEINMANN EXTENSION PIN

with trocar point; triangular shaft

Ø 4.0mm; 14.0cm (5.6")



SKU: 212-348-14

Categories: [Open Surgery](#), [Orthopaedics](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Steinmann pins
Aesculap	KX114
KLS Martin	22-562-14
AUTHORS	STEINMANN

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Steinmann pins
Aesculap	KX114
KLS Martin	22-562-14
AUTHORS	STEINMANN