

STILLE DRESSING FORCEPS

straight; serrated

15.0cm (6"); width of jaw 2.7 +/-0.2

SKU: 020-008-15
 Categories: [Open Surgery](#), [Forceps](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dressing Forceps , Tissue Forceps
Aesculap	BD067
KLS Martin	12-130-15
AUTHORS	STANDARD
Medicon	06.44.15
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dressing Forceps , Tissue Forceps
Aesculap	BD067
KLS Martin	12-130-15

AUTHORS

[STANDARD](#)

Medicon

06.44.15

SPECIFICATIONS

[serrated](#), [straight](#)