

STORAGE AND STERILISATION RACK

for 10 laparoscopy instruments Ø 3.0mm
 Ø 5.0mm or Ø 10.0mm; 525x245x175mm
 incl. 2 silicon holders; for container 580x280x200mm



SKU: 499-101-20
Categories: [Endoscopy](#), [Laparoscopy](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)
RANGE [Accessories](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)
RANGE [Accessories](#)