

## SUBRAMANIAN PEDIATRIC VASCULAR CLAMP

left angled; graduated; atraumatic

16.0cm (6.2")




**SKU:** 142-631-02

**Categories:** [Open Surgery](#), [Thorax](#), [Cardiovascular Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Vascular Clamps</a>
<b>KLS Martin</b>	24-563-16
<b>AUTHORS</b>	<a href="#">SUBRAMANIAN</a>
<b>SPECIFICATIONS</b>	<a href="#">angled</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Vascular Clamps</a>
<b>KLS Martin</b>	24-563-16
<b>AUTHORS</b>	<a href="#">SUBRAMANIAN</a>
<b>SPECIFICATIONS</b>	<a href="#">angled</a>