

SURGICAL SUCTION TUBE

tapered; Ø 4.0mm; 18.0cm

SKU: 085-019-40

Categories: [Open Surgery](#), [Puncture / Special Needles](#)



ZUSÄTZLICHE INFORMATIONEN

| | |
|-------------------|-------------------------------|
| UNIT | pcs. |
| TYPE | Suction Tubes |
| KLS Martin | 18-504-18 |

ZUSÄTZLICHE INFORMATIONEN

| | |
|-------------------|-------------------------------|
| UNIT | pcs. |
| TYPE | Suction Tubes |
| KLS Martin | 18-504-18 |