

## SURGICAL SUCTION TUBE

tapered; Ø 6.0mm; 18.0cm

**SKU:** 085-019-60

**Categories:** [Open Surgery](#), [Puncture / Special Needles](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Suction Tubes</a>
<b>KLS Martin</b>	18-506-18

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Suction Tubes</a>
<b>KLS Martin</b>	18-506-18