

TAYLOR DRESSING FORCEPS

bayonet shaped; serrated

18.5cm (7"); width of jaw 1.5 +/-0.2

SKU: 020-117-18

Categories: [Open Surgery](#), [Forceps](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dressing Forceps , Tissue Forceps
KLS Martin	12-183-17
AUTHORS	BONNEY
Medicon	07.26.18
SPECIFICATIONS	bajonett , serrated



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dressing Forceps , Tissue Forceps
KLS Martin	12-183-17
AUTHORS	BONNEY
Medicon	07.26.18

SPECIFICATIONS

[bajonett, serrated](#)