

TC ADSON DRESSING FORCEPS

straight; serrated; jaw 1,5+/-0.2

15.0cm (6")

SKU: 021-070-15
Categories: [Open Surgery](#), [Forceps](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dressing Forceps , Tissue Forceps
MATERIAL	Tungsten Carbide(TC)
KLS Martin	12-250-15
AUTHORS	ADSON , MICRO
Medicon	06.20.65
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dressing Forceps , Tissue Forceps
MATERIAL	Tungsten Carbide(TC)
KLS Martin	12-250-15

AUTHORS

[ADSON](#), [MICRO](#)

Medicon

06.20.65

SPECIFICATIONS

[serrated](#), [straight](#)