


TC BABY-CRILE-WOOD NEEDLE HOLDER

straight; serrated

15.0cm (6")



SKU: 051-062-15

Categories: [Open Surgery](#), [Suture](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM013
MATERIAL	Tungsten Carbide(TC)
AUTHORS	CRILE , WOOD
Medicon	10.20.64
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM013
MATERIAL	Tungsten Carbide(TC)

AUTHORS

[CRILE](#), [WOOD](#)

Medicon

10.20.64

SPECIFICATIONS

[serrated](#), [straight](#)