

TC BARRAQUER MICRO NEEDLE HOLDER

curved; serrated; with catch

18.0cm (7.2")

SKU: 098-035-18

Categories: [Open Surgery](#), [Suture](#)



ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|--------------------------------------|
| UNIT | pcs. |
| TYPE | Micro-Needle Holders |
| MATERIAL | Tungsten Carbide(TC) |
| AUTHORS | BARRAQUER |
| SPECIFICATIONS | curved, with catch |



ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|--------------------------------------|
| UNIT | pcs. |
| TYPE | Micro-Needle Holders |
| MATERIAL | Tungsten Carbide(TC) |
| AUTHORS | BARRAQUER |
| SPECIFICATIONS | curved, with catch |