

TC BARRAQUER MICRO NEEDLE HOLDER

straight; smooth; with catch

14.0cm (5.6")

SKU: 098-030-14

Categories: [Open Surgery](#), [Suture](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Micro-Needle Holders
MATERIAL	Tungsten Carbide(TC)
AUTHORS	BARRAQUER
SPECIFICATIONS	straight , with catch



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Micro-Needle Holders
MATERIAL	Tungsten Carbide(TC)
AUTHORS	BARRAQUER
SPECIFICATIONS	straight , with catch