

TC BARRY STERNAL WIRE TWISTER

straight; serrated

15.0cm (6")




SKU: 051-240-15

Categories: [Open Surgery](#), [Suture](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Wire twisting Forceps
MATERIAL	Tungsten Carbide(TC)
KLS Martin	22-500-01
AUTHORS	BARRY STERNAL
Medicon	14.40.65
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Wire twisting Forceps
MATERIAL	Tungsten Carbide(TC)
KLS Martin	22-500-01

AUTHORS

[BARRY STERNAL](#)

Medicon

14.40.65

SPECIFICATIONS

[serrated](#), [straight](#)