


TC BOZEMANN NEEDLE HOLDER

s-shaped; straight; serrated

24.0cm (9.6")




SKU: 051-185-24

Categories: [Open Surgery](#), [Suture](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM101
MATERIAL	Tungsten Carbide(TC)
KLS Martin	20-659-24
AUTHORS	BOZEMANN
Medicon	10.37.74
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM101

MATERIAL	Tungsten Carbide(TC)
KLS Martin	20-659-24
AUTHORS	BOZEMANN
Medicon	10.37.74
SPECIFICATIONS	serrated , straight