

TC BOZEMANN NEEDLE HOLDER

s-shaped; straight; serrated

26.0cm (10.4")



SKU: 051-185-26

Categories: [Open Surgery](#), [Suture](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
MATERIAL	Tungsten Carbide(TC)
KLS Martin	20-659-26
AUTHORS	BOZEMANN
Medicon	10.37.76
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
MATERIAL	Tungsten Carbide(TC)
KLS Martin	20-659-26

AUTHORS

[BOZEMANN](#)

Medicon

10.37.76

SPECIFICATIONS

[serrated](#), [straight](#)