

## TC BOZEMANN NEEDLE HOLDER

s-shaped; straight; serrated

20.0cm (8")



**SKU:** 051-185-20

**Categories:** [Open Surgery](#), [Suture](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Needle Holders</a>
<b>Aesculap</b>	BM100
<b>MATERIAL</b>	<a href="#">Tungsten Carbide(TC)</a>
<b>KLS Martin</b>	20-659-20
<b>AUTHORS</b>	<a href="#">BOZEMANN</a>
<b>Medicon</b>	10.37.70
<b>SPECIFICATIONS</b>	<a href="#">serrated</a> , <a href="#">straight</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Needle Holders</a>
<b>Aesculap</b>	BM100

<b>MATERIAL</b>	<a href="#">Tungsten Carbide(TC)</a>
<b>KLS Martin</b>	20-659-20
<b>AUTHORS</b>	<a href="#">BOZEMANN</a>
<b>Medicon</b>	10.37.70
<b>SPECIFICATIONS</b>	<a href="#">serrated</a> , <a href="#">straight</a>