

TC CORWIN WIRE TWISTER

straight; serrated
 16.0cm (6.4")



SKU: 051-236-01
Categories: [Open Surgery](#), [Suture](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Wire twisting Forceps
MATERIAL	Tungsten Carbide(TC)
KLS Martin	22-502-01
AUTHORS	CORWIN
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Wire twisting Forceps
MATERIAL	Tungsten Carbide(TC)
KLS Martin	22-502-01
AUTHORS	CORWIN

SPECIFICATIONS

[serrated](#), [straight](#)