


## TC CRILE-WOOD NEEDLE HOLDER

straight; serrated

18.0cm (7.2")



**SKU:** 051-060-18

**Categories:** [Open Surgery](#), [Suture](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Needle Holders</a>
<b>Aesculap</b>	BM017
<b>MATERIAL</b>	<a href="#">Tungsten Carbide(TC)</a>
<b>KLS Martin</b>	20-646-18
<b>AUTHORS</b>	<a href="#">CRILE</a> , <a href="#">WOOD</a>
<b>Medicon</b>	10.20.68
<b>SPECIFICATIONS</b>	<a href="#">serrated</a> , <a href="#">straight</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Needle Holders</a>
<b>Aesculap</b>	BM017

<b>MATERIAL</b>	<a href="#">Tungsten Carbide(TC)</a>
<b>KLS Martin</b>	20-646-18
<b>AUTHORS</b>	<a href="#">CRILE</a> , <a href="#">WOOD</a>
<b>Medicon</b>	10.20.68
<b>SPECIFICATIONS</b>	<a href="#">serrated</a> , <a href="#">straight</a>