


TC CRILE-WOOD NEEDLE HOLDER

straight; serrated

20.0cm (8")



SKU: 051-060-20

Categories: [Open Surgery](#), [Suture](#)

ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|---|
| UNIT | pcs. |
| TYPE | Needle Holders |
| Aesculap | BM019 |
| MATERIAL | Tungsten Carbide(TC) |
| KLS Martin | 20-638-20 |
| AUTHORS | CRILE , WOOD |
| Medicon | 10.20.70 |
| SPECIFICATIONS | serrated , straight |



ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------|--------------------------------|
| UNIT | pcs. |
| TYPE | Needle Holders |
| Aesculap | BM019 |

| | |
|-----------------------|---|
| MATERIAL | Tungsten Carbide(TC) |
| KLS Martin | 20-638-20 |
| AUTHORS | CRILE , WOOD |
| Medicon | 10.20.70 |
| SPECIFICATIONS | serrated , straight |