

TC CUSHING DRESSING FORCEPS

straight; serrated

17.5cm (7")

SKU: 021-092-17
Categories: [Open Surgery](#), [Forceps](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dressing Forceps , Tissue Forceps
Aesculap	BD160
MATERIAL	Tungsten Carbide(TC)
KLS Martin	12-254-17
AUTHORS	CUSHING
Medicon	06.60.68
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dressing Forceps , Tissue Forceps
Aesculap	BD160

MATERIAL	Tungsten Carbide(TC)
KLS Martin	12-254-17
AUTHORS	CUSHING
Medicon	06.60.68
SPECIFICATIONS	serrated , straight