

TC CUSHING-TAYLOR DRESSING FORCEPS

bayonet shaped; serrated

18.0cm (7.2")

SKU: 021-117-18

Categories: [Open Surgery](#), [Forceps](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dressing Forceps , Tissue Forceps
MATERIAL	Tungsten Carbide(TC)
KLS Martin	12-256-18
AUTHORS	WANGENSTEEN
Medicon	07.26.68
SPECIFICATIONS	bajonett , serrated



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dressing Forceps , Tissue Forceps
MATERIAL	Tungsten Carbide(TC)
KLS Martin	12-256-18

AUTHORS

[WANGENSTEEN](#)

Medicon

07.26.68

SPECIFICATIONS

[bajonett](#), [serrated](#)