


TC DE BAKEY NEEDLE HOLDER

straight; serrated

20.0cm (8")



SKU: 051-114-20

Categories: [Open Surgery](#), [Suture](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM034
MATERIAL	Tungsten Carbide(TC)
KLS Martin	20-646-20
AUTHORS	DE BAKEY
Medicon	10.20.70
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM034

MATERIAL	Tungsten Carbide(TC)
KLS Martin	20-646-20
AUTHORS	DE BAKEY
Medicon	10.20.70
SPECIFICATIONS	serrated , straight