


TC DE BAKEY NEEDLE HOLDER

straight; serrated

31.0cm (12.4")



SKU: 051-114-31

Categories: [Open Surgery](#), [Suture](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM038
MATERIAL	Tungsten Carbide(TC)
KLS Martin	20-646-31
AUTHORS	DE BAKEY
Medicon	10.20.80
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM038

MATERIAL	Tungsten Carbide(TC)
KLS Martin	20-646-31
AUTHORS	DE BAKEY
Medicon	10.20.80
SPECIFICATIONS	serrated , straight