

## TC DERF NEEDLE HOLDER

straight; serrated

12.0cm (4.8")



**SKU:** 051-010-12

**Categories:** [Open Surgery](#), [Suture](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Needle Holders</a>
<b>Aesculap</b>	BM008
<b>MATERIAL</b>	<a href="#">Tungsten Carbide(TC)</a>
<b>KLS Martin</b>	20-632-12
<b>AUTHORS</b>	<a href="#">DERF</a>
<b>Medicon</b>	10.12.62
<b>SPECIFICATIONS</b>	<a href="#">serrated</a> , <a href="#">straight</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Needle Holders</a>
<b>Aesculap</b>	BM008

<b>MATERIAL</b>	<a href="#">Tungsten Carbide(TC)</a>
<b>KLS Martin</b>	20-632-12
<b>AUTHORS</b>	<a href="#">DERF</a>
<b>Medicon</b>	10.12.62
<b>SPECIFICATIONS</b>	<a href="#">serrated</a> , <a href="#">straight</a>