

TC FINOCHIETTO NEEDLE HOLDER

angled; serrated

27.0cm (10.8")




SKU: 051-211-27

Categories: [Open Surgery](#), [Suture](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM092
MATERIAL	Tungsten Carbide(TC)
KLS Martin	20-653-27
AUTHORS	FINOCHIETTO
Medicon	10.43.77
SPECIFICATIONS	angled , serrated



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM092

MATERIAL	Tungsten Carbide(TC)
KLS Martin	20-653-27
AUTHORS	FINOCHIETTO
Medicon	10.43.77
SPECIFICATIONS	angled , serrated