

## TC FINOCHIETTO NEEDLE HOLDER

angled; serrated

27.0cm (10.8")




**SKU:** 051-211-27

**Categories:** [Open Surgery](#), [Suture](#)

### ZUSÄTZLICHE INFORMATIONEN

|                       |   |
|-----------------------|---|
| <b>UNIT</b>           | <a href="#">pcs.</a>                              |
| <b>TYPE</b>           | <a href="#">Needle Holders</a>                    |
| <b>Aesculap</b>       | BM092   |
| <b>MATERIAL</b>       | <a href="#">Tungsten Carbide(TC)</a>              |
| <b>KLS Martin</b>     | 20-653-27   |
| <b>AUTHORS</b>        | <a href="#">FINOCHIETTO</a>                       |
| <b>Medicon</b>        | 10.43.77  |
| <b>SPECIFICATIONS</b> | <a href="#">angled</a> , <a href="#">serrated</a> |



### ZUSÄTZLICHE INFORMATIONEN

|                 |                                |
|-----------------|--------------------------------|
| <b>UNIT</b>     | <a href="#">pcs.</a>           |
| <b>TYPE</b>     | <a href="#">Needle Holders</a> |
| <b>Aesculap</b> | BM092                          |

|                       |   |
|-----------------------|---|
| <b>MATERIAL</b>       | <a href="#">Tungsten Carbide(TC)</a>              |
| <b>KLS Martin</b>     | 20-653-27   |
| <b>AUTHORS</b>        | <a href="#">FINOCHIETTO</a>                       |
| <b>Medicon</b>        | 10.43.77  |
| <b>SPECIFICATIONS</b> | <a href="#">angled</a> , <a href="#">serrated</a> |