

TC FLAT NOSE PLIER DOUBLE CONVERTIBLE

fig. 2; parallel grip

18.0cm (7.2")


SKU: 190-741-18

Categories: [Open Surgery](#), [Orthopaedics](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Flat Noise Pliers
Aesculap	LX182
MATERIAL	Tungsten Carbide(TC)
KLS Martin	22-511-02



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Flat Noise Pliers
Aesculap	LX182
MATERIAL	Tungsten Carbide(TC)
KLS Martin	22-511-02