

TC FLAT NOSE PLIER DOUBLE CONVERTIBLE

fig. 3; parallel grip

18.0cm (7.2")



SKU: 190-742-18

Categories: [Open Surgery](#), [Orthopaedics](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Flat Noise Pliers
Aesculap	LX183
MATERIAL	Tungsten Carbide(TC)
KLS Martin	22-511-03



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Flat Noise Pliers
Aesculap	LX183
MATERIAL	Tungsten Carbide(TC)
KLS Martin	22-511-03