

TC GERALD DRESSING FORCEPS

straight; serrated

18.0cm (7.2")

SKU: 021-102-18
 Categories: [Open Surgery](#), [Forceps](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dressing Forceps , Tissue Forceps
MATERIAL	Tungsten Carbide(TC)
AUTHORS	CUSHING , TAYLOR
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dressing Forceps , Tissue Forceps
MATERIAL	Tungsten Carbide(TC)
AUTHORS	CUSHING , TAYLOR
SPECIFICATIONS	serrated , straight