

TC HALSEY NEEDLE HOLDER

straight; smooth, for left hand use

13.0cm (5.2")

SKU: 051-030-13/LH

Categories: [Open Surgery](#), [Suture](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
MATERIAL	Tungsten Carbide(TC)
AUTHORS	HALSEY
SPECIFICATIONS	straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
MATERIAL	Tungsten Carbide(TC)
AUTHORS	HALSEY
SPECIFICATIONS	straight