

TC HALSEY NEEDLE HOLDER

straight; serrated

13.0cm (5.2")



SKU: 051-040-13

Categories: [Open Surgery](#), [Suture](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM012
MATERIAL	Tungsten Carbide(TC)
KLS Martin	20-637-13
AUTHORS	HALSEY
Medicon	10.12.63
SPECIFICATIONS	serrated , straight

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM012

MATERIAL	Tungsten Carbide(TC)
KLS Martin	20-637-13
AUTHORS	HALSEY
Medicon	10.12.63
SPECIFICATIONS	serrated , straight