

TC HALSEY NEEDLE HOLDER

straight; serrated, for left hand use

13.0cm (5.2")



SKU: 051-040-13/LH

Categories: [Open Surgery](#), [Suture](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
MATERIAL	Tungsten Carbide(TC)
AUTHORS	HALSEY
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
MATERIAL	Tungsten Carbide(TC)
AUTHORS	HALSEY
SPECIFICATIONS	serrated , straight