

TC HALSEY NEEDLE HOLDER

straight; smooth

13.0cm (5.2")




SKU: 051-030-13

Categories: [Open Surgery](#), [Suture](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM010
MATERIAL	Tungsten Carbide(TC)
KLS Martin	20-636-13
AUTHORS	HALSEY
Medicon	10.10.63
SPECIFICATIONS	straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM010

MATERIAL	Tungsten Carbide(TC)
KLS Martin	20-636-13
AUTHORS	HALSEY
Medicon	10.10.63
SPECIFICATIONS	straight