

## TC HALSEY NEEDLE HOLDER

straight; smooth

13.0cm (5.2")




**SKU:** 051-030-13

**Categories:** [Open Surgery](#), [Suture](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Needle Holders</a>
<b>Aesculap</b>	BM010
<b>MATERIAL</b>	<a href="#">Tungsten Carbide(TC)</a>
<b>KLS Martin</b>	20-636-13
<b>AUTHORS</b>	<a href="#">HALSEY</a>
<b>Medicon</b>	10.10.63
<b>SPECIFICATIONS</b>	<a href="#">straight</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Needle Holders</a>
<b>Aesculap</b>	BM010

<b>MATERIAL</b>	<a href="#">Tungsten Carbide(TC)</a>
<b>KLS Martin</b>	20-636-13
<b>AUTHORS</b>	<a href="#">HALSEY</a>
<b>Medicon</b>	10.10.63
<b>SPECIFICATIONS</b>	<a href="#">straight</a>